



Reference: PTXBC024021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MBD2; dipeptidase 2		
For research use only.			

GENE DESCRIPTION	dipeptidase 2
GENE ID	64174
ORIGIN SPECIES	Human
ACCESSIONS	BC024021
VECTOR	pDONR223
NCBI SYMBOL	DPEP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]