

SIX2-SIX homeobox 2 Gene

Reference: PTXBC024033

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	homeobox protein SIX2; sine oculis homeobox homolog 2; SIX homeobox 2		
For research use only.			

GENE DESCRIPTION	SIX homeobox 2
GENE ID	10736
ORIGIN SPECIES	Human
ACCESSIONS	BC024033
VECTOR	pDONR223
NCBI SYMBOL	SIX2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtccatgctccaccttcggcttcacgcaggagcaagtggcgtgctgtgcagcagggcgcaacatcgagcggctggcgcttcctgtgtgctgctgcccgcctgcgagcacttcacaagaatgaaagcgtgctcaaggccaaggccgtggtggccttcacc
gcggcaacttcgcgagctctacaagatcctggagagccaccagtctcgcgcacacccacgcaagctgcagcagctgtggtcgaagcacaactacatcgaggcggagaagctgcgcggccgaccctggcgccgtgggcaaataccgctgctgcgcaattcccgctgcc
gcgctccatctgggacggcgaggagaccagctactgtctcaaggaaaagatcgagcgtgctgctgcgagtggttacgcgcacacccctacccctcaccgcgagaagcgtgagctgacggagccacgggctcaccaccacacaggtcagcaactgttcaagaaccggcgg
cagcgcgaccggcgccgagccaaggaaagggaacaacgagaactccaattctaacagccacaacccctgaatggcagcggcgaagtcgggtgttaggcagctcggaggatgagaagactccatcggggacgccaagaccactcatccagccccgactgctcctcagc
ccggcggccctgggtgctgctcctgcacagcctggggccaccctccggggcccgagcgcagtgccagtgccgggtgccagggcggaggtggagcggaccactgcaacaccacatggcctgcaggactccatcctcaacccatgtcagccaacctgtgacctgggctcctag