



Reference: PTXBC024039

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AFP; RNF95; tripartite motif-containing protein 26; RING finger protein 95; acid finger protein; widely expressed acid zinc finger protein; zinc finger protein 173; tripartite motif containing 26		

For research use only.

GENE DESCRIPTION	tripartite motif-containing 26
GENE ID	7726
ORIGIN SPECIES	Human
ACCESSIONS	BC024039
VECTOR	pDONR223
NCBI SYMBOL	TRIM26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr