

POLR2J-polymerase (RNA) II (DNA directed) polypeptide J, 13.3kDa Gene

Reference: PTXBC024165

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	POLR2J1; RPB11; RPB11A; RPB11m; hRPB14; DNA-directed RNA polymerase II subunit RPB11-a; DNA-directed RNA polymerase II subunit J-1; RNA polymerase II 13.3 kDa subunit; RNA polymerase II subunit B11-a; polymerase (RNA) II (DNA directed) polypeptide J, 13.3		

For research use only.

GENE DESCRIPTION	polymerase (RNA) II (DNA directed) polypeptide J, 13.3kDa
GENE ID	5439
ORIGIN SPECIES	Human
ACCESSIONS	BC024165
VECTOR	pDONR223
NCBI SYMBOL	POLR2J
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaacgcccctccagccttcgagtcgttctgtcttcgagggcgagaagaagatcaccattaacaaggacaccaaggtacccaatgctgtttattcaccatcaacaagaagaccacacactgggaacatcattaatacacaactcctaaagaccgcaagtgcatttgctggctac
aaagtcccccacccttgagcacaagaatcatccagtgagaccacgcccgaactacgccccaggaagcctttaccaacgccatcaccgacccatcatcagtgagctgtccctgctggaggagcgtttcgggtggccataaaagacaagcaggaaggaattgagtag