

C1orf150-chromosome 1 open reading frame 150 Gene

Reference: PTXBC024174

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C1orf150; germinal center-associated signaling and motility-like protein; germinal center associated signaling and motility like		

For research use only.

GENE DESCRIPTION	chromosome 1 open reading frame 150
GENE ID	148823
ORIGIN SPECIES	Human
ACCESSIONS	BC024174
VECTOR	pDONR223
NCBI SYMBOL	C1orf150
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaattatctcctgcgaaaactcagttgcctgggagagaatcaaaagaagcccaagaaaggaaacccagatgaggaaagaaaacggcaggaaatgactacatttgaagaaaactcaagatcaagataagaaaagccaagaaagtttcacacttctaatcaggaaaacgagaat
ggcagtggttctgaagaagtggtctacactgtcattaacacatccccatcagagatcctccctgagctccaatgatgatgctatgagaacattgactccctcacaagaaagtgagacagtttagagaaggtcagagacagaatatgccctcttaggacttctgttagtaggccttgttcc
tgcacccatgagcatgattatgaagttgtgttccacactaa