

ROBLD3-roadblock domain containing 3 Gene

Reference: PTXBC024190

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ROBLD3; ENDAP; HSPC003; MAPBPIP; MAPKSP1AP; Ragulator2; p14; ragulator complex protein LAMTOR2; MAPBP-interacting protein; MAPKSP1 adaptor protein; endosomal adaptor protein p14; late endosomal/lysosomal Mp1-interacting protein; late endosomal/lysosomal a		

For research use only.

GENE DESCRIPTION	roadblock domain containing 3
GENE ID	28956
ORIGIN SPECIES	Human
ACCESSIONS	BC024190
VECTOR	pDONR223
NCBI SYMBOL	ROBLD3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgcgccccaggcttggaccaggtgctaagccaagccaacactggaggcgtccagagaccctgtctgaataacgaggatcactgctggcctactctggttacgggacactgacgccccgggtcaccgctgccatagccagtaacatctggccgctacgaccggaacg
ggaaccaagcgtttaatgaagacaatctcaattcatcctcatggactgcatggaggccgtgtagccatcaccgagtgccaaccttctgtgtgtatgtatgccaaggagaccgtggcctttggaatgctcaaggccaaggccagccttgggtgactgacctggagagccctcac
ccaagtggcgcgcatctaa