

SUPT5H-suppressor of Ty 5 homolog (*S. cerevisiae*) Gene

Reference: PTXBC024203

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SPT5; SPT5H; Tat-CT1; transcription elongation factor SPT5; DRB sensitivity-inducing factor 160 kDa subunit; DRB sensitivity-inducing factor large subunit; DSIF large subunit; DSIF p160; Tat-cotransactivator 1 protein; hSPT5; suppressor of Ty 5 homolog; S		

For research use only.

GENE DESCRIPTION	suppressor of Ty 5 homolog (S. cerevisiae)
GENE ID	6829
ORIGIN SPECIES	Human
ACCESSIONS	BC024203
VECTOR	pDONR223
NCBI SYMBOL	SUPT5H
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr