

M6PR-mannose-6-phosphate receptor (cation dependent) Gene

Reference: PTXBC024206

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD-M6PR; CD-MPR; MPR 46; MPR-46; MPR46; SMPR; cation-dependent mannose-6-phosphate receptor; 46-kDa mannose 6-phosphate receptor; CD Man-6-P receptor; Mr 46,000 Man6PR; small mannose 6-phosphate receptor; mannose-6-phosphate receptor, cation dependent		

For research use only.

GENE DESCRIPTION	mannose-6-phosphate receptor (cation dependent)
GENE ID	4074
ORIGIN SPECIES	Human
ACCESSIONS	BC024206
VECTOR	pDONR223
NCBI SYMBOL	M6PR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttcccttctacagctgctggaggactggactgctactactactcctggctgtggcagtgagagaatcctggcagacagaagaaaaacttgcgacttggtaggagaaaaagggtaagagtcagagaaagagttggctctagtgaaaggctgaaccactgtttaataaaagctttga
gagcactgtggccagggttcagacacatacatcttcagggtgtgccgggaagctggcaaccacactctggggcaggcctgtgcaaatcaacaaagtaafgggaaggagacagtggttagggagactcaacgagactcacatcttcaacggaagtaattggatcatgctgat
ctataaaggggtgatgaatatgacaaccactgtggcaaggagcagcgtcgtgagtgatgatctcctgcaatcgacacaccctagcggacaatttaacccctgtgtctgaggagcgtggcaagtcgaagattgtttctaccttttgatggatagcagcctggcctgttcaccaga
gatctccacacagtggtgggttcacatttgcacgtttgcacactggtgtgtttatgtttgttgggggttcctataccagcagctgttagtgggagccaaaggaaatggagcagtttcccacttagccttctggcaggatctggcaacctggtagcagatggcgtgactttgtctgcc
gttctaaacctgaaatgtgcctgcagcatatcgtgtgtgtggggtgaccagctggggaggagtcagaagaagggtgaccatttattaccaatgtag