

LOC202459-similar to RIKEN cDNA 2310008M10 Gene

Reference: PTXBC024224

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	similar to RIKEN cDNA 2310008M10
GENE ID	202459
ORIGIN SPECIES	Human
ACCESSIONS	BC024224
VECTOR	pDONR223
NCBI SYMBOL	LOC202459
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgattgtatgctctgatgggtgtcttactctctcactggaggaataattatgatgtattgtgaacctccaagcattggctctatgactgatgaacacgggcacagaggccagtagctttcttgcccaacagagtaaatgaacaatgtattatggaaggactgcacccagcttctgt
ttacaataggaggttaggtttcattctctggaccgatggatgcaccaaatatcccaaaactcaatagattctcttattcattggattcgtttgtgctctattgagctttttcatggctagagtattcatgagaatgaaactgccgggctatctgatgggttag