



Reference: PTXBC024229

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EP3-I; EP3-II; EP3-III; EP3-IV; EP3-VI; EP3e; PGE2-R; prostaglandin E2 receptor EP3 subtype; PGE receptor, EP3 subtype; PGE2 receptor EP3 subtype; prostaglandin E receptor EP3 subtype 3 isoform; prostaglandin E receptor 3 (subtype EP3); prostaglandin rece		

GENE DESCRIPTION	prostaglandin E receptor 3 (subtype EP3)
GENE ID	5733
ORIGIN SPECIES	Human
ACCESSIONS	BC024229
VECTOR	pDONR223
NCBI SYMBOL	PTGER3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

Page 1/1