

MGC27345-hypothetical protein MGC27345 Gene

Reference: PTXBC024231

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein MGC27345		

For research use only.

GENE DESCRIPTION	hypothetical protein MGC27345
GENE ID	157247
ORIGIN SPECIES	Human
ACCESSIONS	BC024231
VECTOR	pDONR223
NCBI SYMBOL	MGC27345
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccacaaagcgaagaacacctcttggaggcagagccctcagccctatccaggccagcagtgccacttgaaatctgaaagcagccctccaccaaaccagcaagatggaaaacctccctgggttggcaagtggccctgggaagtacatgccccatctctgcaggatagcagtggaaggttctgcccagtgccagctccaccaccagatgggtactcaagactccaagcctgtcaatgggtggaaacttctggcatttctgtcagaagggtgagtttag