

SCGB3A2-secretoglobin, family 3A, member 2 Gene

Reference: PTXBC024232

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LU103; PNSP1; UGRP1; pnSP-1; secretoglobin family 3A member 2; pneumo secretory protein 1; uteroglobin-related protein 1		

For research use only.

GENE DESCRIPTION	secretoglobin, family 3A, member 2
GENE ID	117156
ORIGIN SPECIES	Human
ACCESSIONS	BC024232
VECTOR	pDONR223
NCBI SYMBOL	SCGB3A2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagctgtaactatctctgctggtgacctcagcctttgtagttactctgctactgacctctcatcaacaagtgcccttctgttgacaagttggcactttacctctggacaacattctccctttatggatccattaaagcttctctgaaaactctgggcatttctgttgacacctgtgg
aggggctaaggaaagtgtgtaaatgagctgggaccagaggcttctgaagctgtgaagaaactgctggagcgctatcacacttggtgtga