

SALL2-sal-like 2 (Drosophila) Gene

Reference: PTXBC024245

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein SALL2; COLB; HSAL2; Sal-2; ZNF795; p150(Sal2); sal-like protein 2; zinc finger protein 795; zinc finger protein Spalt-2; spalt like transcription factor 2		

For research use only.

GENE DESCRIPTION	sal-like 2 (Drosophila)
GENE ID	6297
ORIGIN SPECIES	Human
ACCESSIONS	BC024245
VECTOR	pDONR223
NCBI SYMBOL	SALL2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgacgaatccgagaggagctctcgtcgggtgcccgcggggagccggcagagctcggaggtgatctagcaggaggatcaccccaagctgtgccaaagtgtgcgcacaattcactgacccaactgaattctcgtcccaccagaacgcatttactgacctcctgt
aatggatgaattggggccaggagaaacccaacaactctcggcctcctctgaacccggcctgaggggtcacaataatctcagggtcatggacacagagcatagcaacccccagattctgggtcctcctgtgccacggatcccacctggggccagagaggagaggagagct
ccaggacattctcgtgctgccacagaaccagatgtggcattcctgtcaaatggcctgcccataagccctggaaftccagctccacctccactaccatccaagcctggccccaccagtgtgttggcctaggaaactgtggtgggaaggtgcctccaacaatgggatccaggaa
gccaaggagaagacagccccctctatttcagcctcctgcacccaaggcagtgctga