

RAD54L2-RAD54-like 2 (*S. cerevisiae*) Gene

Reference: PTXBC024298

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ARIP4; HSPC325; SRISNF2L; helicase ARIP4; AR interacting protein 4; RAD54-like protein 2; androgen receptor-interacting protein 4; RAD54-like 2 (S. cerevisiae)		

For research use only.

GENE DESCRIPTION	RAD54-like 2 (S. cerevisiae)
GENE ID	23132
ORIGIN SPECIES	Human
ACCESSIONS	BC024298
VECTOR	pDONR223
NCBI SYMBOL	RAD54L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1