

PHPT1-phosphohistidine phosphatase 1 Gene

Reference: PTXBC024648

| | | | |
|----------------|--|---------------------|----------------|
| Product type: | DNA | Storage Condition: | -20°C |
| Size/Quantity: | 2 ug | Delivery Condition: | +4°C, Blue Ice |
| Synonyms | CGI-202; HEL-S-132P; HSPC141; PHP; PHP14; 14 kDa phosphohistidine phosphatase; epididymis secretory sperm binding protein Li 132P; protein histidine phosphatase; protein janus-A homolog; sex-regulated protein janus-a; phosphohistidine phosphatase 1 | | |

For research use only.

| | |
|-------------------------|---|
| GENE DESCRIPTION | phosphohistidine phosphatase 1 |
| GENE ID | 29085 |
| ORIGIN SPECIES | Human |
| ACCESSIONS | BC024648 |
| VECTOR | pDONR223 |
| NCBI SYMBOL | PHPT1 |
| LENGTH | ProteoGenix bp |
| SEQUENCE PRIMERS | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

ORF (Only for reference, product provided without STOP codon)

atggcgggtggcggacctcgtctcattctgatgtggacatcgactccgacggcgcttcaagtatgtgctgatccgagtcactcggctccccgctccgggctccggctgcagagagcaaggagatcgtgcggctacaagtggctgagtaccatcgggacatctacgacaaaatg
tcggcgacatgcagaagcaaggctgcgactgtgagtgctggcgggcgcatctccaccagatcaggacaaagaattcacgtgtacggctattcattgcctatggtcctgcccagcagccatttcaactgagaaaaatcaagcaagtaaccccgactacgaggtcacctg
ggctaacgacggctactga