

PTPN11-protein tyrosine phosphatase, non-receptor type 11 Gene

Reference: PTXBC025181

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BPTP3; CFC; JMML; METCDS; NS1; PTP-1D; PTP2C; SH-PTP2; SH-PTP3; SHP2; tyrosine-protein phosphatase non-receptor type 11; PTP-2C; protein-tyrosine phosphatase 1D; protein-tyrosine phosphatase 2C; protein tyrosine phosphatase, non-receptor type 11		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, non-receptor type 11
GENE ID	5781
ORIGIN SPECIES	Human
ACCESSIONS	BC025181
VECTOR	pDONR223
NCBI SYMBOL	PTPN11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcttggtataaaggatgatcgaggttctcagacgagtttactttacatgagatggaatcaggcagagaggctggatgatggagaagctcaggtgaagtgtttataaaaaaagggtgtggaaggaaagtccaagaggtggtttctga