

AIG1-androgen-induced 1 Gene

Reference: PTXBC025278

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AIG-1; dJ95L4.1; androgen-induced gene 1 protein; androgen induced 1		
For research use only.			

GENE DESCRIPTION	androgen-induced 1
GENE ID	51390
ORIGIN SPECIES	Human
ACCESSIONS	BC025278
VECTOR	pDONR223
NCBI SYMBOL	AIG1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgcttgcctcccgaggtgctgcggatggcaatcctgctgtcttactgctctatcctgtgtgaactacaaggccatcgaaatgcctcacaccagacctacggaggagctggaaattcctgacgttcattgatctggttatccaggctgtctttttggcatctgtgctgaactgtcttcc
agctcttctgactcgaggagtgaggaccaggagcaagagaggcagctcaagaagctcatctctcctgggactggatgttagctgtgttggccttctctgttgggtttttgtgtagcagtgcttgatcattatgcctatgacagagagatgatataccggaagctgctggataattttatccc
agggtggctgaatcacggaatgcacacgacggttctgcccctttatattatcgagatgaggacatcgcacatcagtatccagcaggagcagcggactaccgcatatgtaccttctctgttggctatatattgggtgtgctgggtgcatcatgtaactggcatgtgggtgtacccttctc
ggaacacattggccaggagccagaatcctcttgggtctacaacctcttaataaacttctgtacgtgctggagaagttctgaacaactatatctgggatacacagaaaagtatggaagaagagaaaagcctaaattggaatga