

RNASE11-ribonuclease, RNase A family, 11 (non-active) Gene

Reference: PTXBC025410

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C14orf6; HEL-S-84p; RAJ1; RNase 11; epididymis secretory sperm binding protein Li 84p; ribonuclease 11; ribonuclease A J1; ribonuclease, RNase A family, 11 (non-active); ribonuclease A family member 11 (inactive)		

For research use only.

GENE DESCRIPTION	ribonuclease, RNase A family, 11 (non-active)
GENE ID	122651
ORIGIN SPECIES	Human
ACCESSIONS	BC025410
VECTOR	pDONR223
NCBI SYMBOL	RNASE11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaccttctctgctgctcagcctggcctgttcttcagaagcatcagaagcacaatgaagataattaagaagaatttacagacgaagagatgcaatgacatggcaaaagtggccaagaaaaacagaccattgagatattaatgaaccgatcctgttagttaaata
ccagcctcagcatgtccaaggatgatgtcttccacattactgacattcagaagttacattataatgacccaagggaacagttcggtaatagacaagagtggtgcaatgacatgacagctctggagaaagtgttcagaagcaaacggatcgtgcaagtggagcaataactccgcag
ctccacagaagtgtatgcgcagggtccacagggcccccagctgcaagttttgacagaatcctggcataagctgctgtgagagcctagaactggaataacagtggtccagttcactacaggcaacaattccccaggtgccaataccatagtggtacctcattagagaagatattgacagtgc
tgacaggtcattctctgatgagctggttagttgtggtcttaagttgtaa