



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

C11orf41-chromosome 11 open reading frame 41 Gene

Reference: PTXBC025668

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0606 protein C11orf41; C11orf41; C11orf69; UPF0606 protein KIAA1549L; Uncharacterized protein C11orf69; KIAA1549 like		

For research use only.

GENE DESCRIPTION	chromosome 11 open reading frame 41
GENE ID	25758
ORIGIN SPECIES	Human
ACCESSIONS	BC025668
VECTOR	pDONR223
NCBI SYMBOL	C11orf41
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccacgcagcagcagcagcagcagctcggaagcagctccccctgactaccacacagagagaagaaggggcctcatcttgcacatcaggctgaggggtgcacacacacttgaactctgctagactggccagagctgccccagccattctcaccaccaccagctattgatacagccgtgggtct
aggaagcgaattctgaacaagcccgccagttgacacaaatctctctcagcagacaaataaacttctgaattgaattgcacatcttactcagctgctactaatctggctcacgccacacacttgcattgataa