

## CCDC33-coiled-coil domain containing 33 Gene

Reference: PTXBC025689

|                |                                                                                                                     |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                 | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | CT61; HP11097; coiled-coil domain-containing protein 33; cancer/testis antigen 61; coiled-coil domain containing 33 |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | coiled-coil domain containing 33                                                                  |
| <b>GENE ID</b>          | 80125                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC025689                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | CCDC33                                                                                            |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgggggagccaggccccccatccccccagcccccaagcgtgccccagccctcattgacagcttctctgagatgaacaactaccggcgccatgcagaagatggcagaggacatcctgtctctcgggagacaggccagcatcctggaaggagagaaccgcatactgaggagccg  
cctggccccagcaggaggaggaagaggggcaggcgaagccagtgaggccccagaacacgggtgtccatgaaagcagaactgctgctgagtgagctggatgaaagaactgaggacagggtgcagcattgcaagaatgagctgattcgaagaatgatcgagagaaggagctgctc  
cttctgtatcaggccccagcagccacagggcgtctgctggaagcagtgaccggcgaagctgcagaagatgaaggcgtgaggagagctgtcggcaccagaagaaggtgatcgagaagatggagcgggtgctggaggacaggctgcaggacaggagcaagccccctcctgaac  
aggcagcagggaagccctacacgggcttccctatgctctcagcctctggccttccctgggttctatgggagagaacctgccggtgaactttactcgtgctgtggcagaaaaacgcgaagctcggacggagctggataagaaccgccaccagcagcccatcattctcagcaaca  
ggccctgcccgtggagccccggggagttgggagcaggaggagacttgacagagaaggctacaagagacgcatggccccagccactcagagtgacagagaccctgccccgacaggatctcctctgtgacttcagacaagttcaacctcctggccaagctggaacacgctcagagcc  
ggatcctgtccctggaagccaggttagaggactcagctcagcgtggggagcagagaagcaggatctggccacacggctgcaggagcaagaaaagggttcaggcaccctcgaactccatcatagaacagccctgtgagtga