



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LPAR2-lysophosphatidic acid receptor 2 Gene

Reference: PTXBC025695

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EDG-4; EDG4; LPA-2; LPA2; lysophosphatidic acid receptor 2; G protein-coupled receptor; LPA receptor 2; LPA receptor EDG4; endothelial differentiation, lysophosphatidic acid G-protein-coupled receptor, 4; lysophosphatidic acid receptor EDG4; lysophosphati		

For research use only.

GENE DESCRIPTION	lysophosphatidic acid receptor 2
GENE ID	9170
ORIGIN SPECIES	Human
ACCESSIONS	BC025695
VECTOR	pDONR223
NCBI SYMBOL	LPAR2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr