



Reference: PTXBC025698

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	acyloxyacyl hydrolase (neutrophil)		
For research use only.			

GENE DESCRIPTION	acyloxyacyl hydrolase (neutrophil)
GENE ID	313
ORIGIN SPECIES	Human
ACCESSIONS	BC025698
VECTOR	pDONR223
NCBI SYMBOL	AOAH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]