

## PLA2G2D-phospholipase A2, group IID Gene

Reference: PTXBC025706

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PLA2IID; SPLASH; sPLA2-IID; sPLA2S; group IID secretory phospholipase A2; GIID sPLA2; phosphatidylcholine 2-acylhydrolase 2D; secretory phospholipase A2s; secretory-type PLA, stroma-associated homolog; phospholipase A2 group IID		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	phospholipase A2, group IID
<b>GENE ID</b>	26279
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC025706
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	PLA2G2D
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atggaactgcaactgctgtgtggctggtgatggctggtgattccaatccagggcggatcctgaacctgaacaagatggtcaagcaagtgactggaaaatgccatcctctactggcctacggctgtcactcggactaggtggcagaggccaacccaagatgccag  
gactggtgctgccagaccatgactgctgctatgaccacctgaagaccaggggtgagcatctacaaggactattacagatacaactttccaggggaacatccactgctcgacaagggaagctggtgtgagcagcagctgtgtcctgtgacaaggaggtggcctctgctggaag  
gcaacctggacacctaccagaagcagactgctgttctactggcggccctgcccggggcagaccctgggtgctag