

TRAV20-T cell receptor alpha variable 20 Gene

Reference: PTXBC025727

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	T cell receptor alpha variable 20
GENE ID	28663
ORIGIN SPECIES	Human
ACCESSIONS	BC025727
VECTOR	pDONR223
NCBI SYMBOL	TRAV20
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaaatgttgagtgatcatagctgtggtcagctgtggtgagtgagaagaccaggtgacgagagtcggagccctgagactccaggaggagagtagcagctcaactgcagttacacagtcagcgggttaagagggtgttctggtataggcaagatc
ctgggaaaggccctgaattctcttcacctgtattcagctggggaagaaaggagaaagaaagctaaaggccacattacaagaaggaagctttctgcacatcacagccctaaacctgaagactcagccattatctgtgctgtcggggaggagggaagccaaggaatctca
tctttggaaaaggcactaaactctctgttaaaccaaatccagaacctgacctgcccgtgtaccagctgagagactctaaatccagtgacaagctgtctgcttaccgattttgattctcaacaaatgtgcacaaagtaaggattctgatgtgtatatcagacaaaactgtgctagac
atgaggtctatggacttcaagagcaacagtgctgtggtcggagcaacaatctgacttgcattgtgcaaacgcttcaacaacagcattatccagaagacaccttctccagcccagaagttcctgtgatgtcaagctggtcgagaaagctttgaaacagatacgaacctaacttca
aaacctgtcagtgattgggtccgaatcctctctgaaagtggccgggttaatctgctcatgacgctgcccgtgtgtccagctga