



Reference: PTXBC025737

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	pyruvate kinase PKLR; PK1; PKL; PKR; PKRL; RPK; R-type/L-type pyruvate kinase; pyruvate kinase 1; pyruvate kinase isozyme R/L; pyruvate kinase isozymes L/R; pyruvate kinase isozymes R/L; pyruvate kinase type L; pyruvate kinase, liver and blood cell; red c		

For research use only.

GENE DESCRIPTION	pyruvate kinase, liver and RBC
GENE ID	5313
ORIGIN SPECIES	Human
ACCESSIONS	BC025737
VECTOR	pDONR223
NCBI SYMBOL	PKLR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]