

ITM2C-integral membrane protein 2C Gene

Reference: PTXBC025742

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BRI3; BRICD2C; E25; E25C; ITM3; integral membrane protein 2C; BRICHOS domain containing 2C; cerebral protein 14; integral membrane protein 3; transmembrane protein BRI3		

For research use only.

GENE DESCRIPTION	integral membrane protein 2C
GENE ID	81618
ORIGIN SPECIES	Human
ACCESSIONS	BC025742
VECTOR	pDONR223
NCBI SYMBOL	ITM2C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaagattagctccagcccgctggcctgcatcaaggcgacaaaggctgacaaggcgtcggcgtcggccctgcggcctcggccacccagagatcctgctgacgcccgttagggaggagcagccccacacatc gatccaagagggggagctcagtgggcggtg
tgctacctgtcgaaggcatggtcgtgctgctcagggcctcgtgttcgctctgtctacatctacagatacttcttctgacagctggcccgagataacttctccgctgtggtgtgctgtatgaggactccctgtcctccaggtccggactcagatggagctggaagagatgaaaaatct
acctcgacgagaactacgagcgcacacgtgctgtgccccagtttggcgcggtgacctgcagacatcatccatgacttccagcggagggggacctacctgcgcagacgtacatcatccaggaggagatggtgtcagggagcatgtcagtgacaaggaggccctgggtcctt
catctaccacctgtgcaacgggaaagacacctaccggctcggcgccgggcaacgcggaggcggatcaacaagcgtggggccaagaactgcaatgccatccgccacttcgagaacaccttctgtgtggagacgctcatctgcgggtgtgtgtga