

FAM19A1-family with sequence similarity 19 (chemokine (C-C motif)-like), member A1 Gene

Reference: PTXBC025746

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM19A1; TAFA-1; TAFA1; chemokine-like protein TAFA-1; family with sequence similarity 19 (chemokine (C-C motif)-like), member A1; family with sequence similarity 19 member A1, C-C motif chemokine like		

For research use only.

GENE DESCRIPTION	family with sequence similarity 19 (chemokine (C-C motif)-like), member A1
GENE ID	407738
ORIGIN SPECIES	Human
ACCESSIONS	BC025746
VECTOR	pDONR223
NCBI SYMBOL	FAM19A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcaatgtctctcgatgtcctgggtcctgtatttgggataagtgcttggcaatgctacttgcctggatcccttcagcacacttccagcagcatcacctgcacagaccagaaggaggacgtgtgaagtatagcagcacaccgatgttgaacaagaatcgattgaggagcgt
cacaacagtaaaagtgttctgtctacttgaaaagtggctggaacaacaagaaccggccttcttgcgtcgtatgcctccatagtgattgggaatgggtgtgagatggagccttgctagaaggagaagaatgaagacactccctgacaattctggatgtgcgcaacaggcaa
caaaattaagaccacgagaattcacccaagaacctaa