

MYOZ1-myozenin 1 Gene

Reference: PTXBC025753

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CS-2; FATZ; MYOZ; myozenin-1; calsarcin-2; filamin-, actinin- and telethonin-binding protein; myozenin 1		

For research use only.

GENE DESCRIPTION	myozenin 1
GENE ID	58529
ORIGIN SPECIES	Human
ACCESSIONS	BC025753
VECTOR	pDONR223
NCBI SYMBOL	MYOZ1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccctctcaggaaacccggcccttaataagaagagaaatccagcaagctgatcatggaactcactggaggtggacaggagagctcaggctgaacctggggcaaaaagatcagtgccaaggatgtgatgttgaggaaactgtcgtcttaccacccgggctccaagatgtt
caaaactgcggcagatgaggggtggagaagttaattatgagaaccaccctgatgtttctctgacagctcaatggatcacttcagaagtctctccaacagtggggggacagctgggcacacagctgggtcagggtatctcatacagcaagagcaacggcagagcggcagccaggcagggg
gcagtggtctctcgggacagatggtctgatcagcagcaccatctgggctctgggtctggagctgggggtacaggtgtcccgccggccagctggcagagggagagctgtggcacagcagggtgtgtgagacaggtcagagaccaggcagcgaggagaagaaaacata
tactgtgttcaagacctatatttcccatgggagcgagccatgggggttgacccccagcaaaaatggaacttggcattgacctgtggcctatggggccaaagctgaacttcccaataataagtcctcaacaggacggcaatgccctatggtgatatgagaaggcctccaaacgatg
acctccagatgcccaagtttgacctggggcccttgctgagtgaacctctgtcctctacaacaaaacctctcaacaggccctcttcaatcgaacctctattccctggctgagctctggggagcctgtagactacaacgtggatattggcatcccttgatggagaacagaggagctgt
ga