

C11orf45-chromosome 11 open reading frame 45 Gene

Reference: PTXBC025756

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	chromosome 11 open reading frame 45		
For research use only.			

GENE DESCRIPTION	chromosome 11 open reading frame 45
GENE ID	219833
ORIGIN SPECIES	Human
ACCESSIONS	BC025756
VECTOR	pDONR223
NCBI SYMBOL	C11orf45
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgacgcggctggctcctcagtcacacctgagtagcagcagctctccgccctggagcagcgtgccatcagctgggagctggacaacgtgctgatgcctagtcccagaatctggccccagtgactccaacaggcaggtctgcctctgtcaggagtgagggtaacacctctcactc
tggaaattctcagctggcaggtatgtcatgccatagtaaccagaacctgtgagctctgtgagctctgccgtctaccccacggctgtggctgtgtgaggtctgccacaacattacctgtcagtcctcaggacaacaaggcagcggccccggcaggagaagagagaactcaatctgca
aggcccatgatagtagagagggccgctgggctacccctcagtgcccatcagcctggttcgggtggtcctaactag