

STOML3-stomatin (EPB72)-like 3 Gene

Reference: PTXBC025760

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Epb7.2l; stomatin-like protein 3; EPB72-like 3; SLP-3; erythrocyte band 7 integral membrane protein; protein 7.2B; stomatin (EPB72)-like 3; stomatin like 3		

For research use only.

GENE DESCRIPTION	stomatin (EPB72)-like 3
GENE ID	161003
ORIGIN SPECIES	Human
ACCESSIONS	BC025760
VECTOR	pDONR223
NCBI SYMBOL	STOML3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atggattcagggtgtctccacctgagaagcaagataaagaatctcgtgggtgcaacaataaacggcctgggtatggtgctgatcctgtttccctctctcctgttggatcattaccctcccatcccatatgatgtgctgaagatcattaaggagatgaacgtgctgttattccgt
ctgggacgatccaagctgacaaaagccaaggccaggtttgatcctggtcctgcatgcatagatgtgtttgcaaaagtgacctccgaacagtactgcaacattcctccacaagagatcctcaccagagactccgtaactcaggtagatggagtgtctattacagaatctatagtgct
gtctcagcagtggttaatgtcaacgatgtccatcaagcaacattctgctgctcaaacactctgagaaatgtcttagggacacagacctgtcccagatcttagctggacgagaagagatcggccatagcatccagacttactgatgatgccaccgaactgtggggatccgggtggcc
cgagtggaaatcaagatgttcgattccctgagctgagagatccatggcagccgagcctgagccaccgggaagcagagccaagctctgagctgaaggagaaatgaatgctccaatccctgaagtcagcctccatggtgctgctgagctcctccatagctctccagctg
cgctacctgagacctgagcaggttagccaccgagaagaattctcagattgtttcctctgcccataataactagaggcattgggtggcgtcagctatgataaccacaagaagctccaataaagcctga
```