

FAM76A-family with sequence similarity 76, member A Gene

Reference: PTXBC025768

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM76A; family with sequence similarity 76 member A		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 76, member A
GENE ID	199870
ORIGIN SPECIES	Human
ACCESSIONS	BC025768
VECTOR	pDONR223
NCBI SYMBOL	FAM76A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggcgctctacgcctgcaccaagtgcaccagcgtcccttcgagggcgctgtctcagggcgagcagctgtgcaaggaatgtcggattgcacacctgttgaagtgcacctactgcaggactgagtagcaggagagagtaaaaccaatacaatatgcaagaaatgtgctcag
aacgtgcagttgtatggaacgccaacacgttcagttatgcaacataattgcagcatttattgggaataaatgccagcgtgcacaaatcagaaaagatggaaccacctattctgtgaacagtgcaggcagcgtgtgcatlgacaggaagatgatagaagaagtagatggg
aaattgctgtgctgctgtgcacactttcatacaaacgggtccttcagaagaccaagagcagagggaacacctgagtagctcttctcgtgctggccaccaggaagagcagtagatgcctgagtggtggtggccattataacagccagaaaacacattctacatctcaattcaaat
gaaatcccaagaaaaagtcgaagtttgagtcacataactaatggagacagcttctcccgacctgctctgactcaccagcagcagcacttgcacatcattgcccaactgaaggagaagtggtgctacctgaagaagatgttgatcaaaaggatcaaatgattttagagaagag
aagaagattacagagttgaagctgattttcagtagcaggaatgcagatgagagccaaaatgaaccagatggagaaaaccacaaaagatcacagaaactgcaggccaaaaccgagagctcctgaagcaggcagctgttctcgaagcaagaagtcagagaagtcagga
gctataacctctcatga