

NPPB-natriuretic peptide precursor B Gene

Reference: PTXBC025785

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	BNP; natriuretic peptides B; brain type natriuretic peptide; gamma-brain natriuretic peptide; natriuretic peptide precursor B; natriuretic protein; natriuretic peptide B		

For research use only.

GENE DESCRIPTION	natriuretic peptide precursor B
GENE ID	4879
ORIGIN SPECIES	Human
ACCESSIONS	BC025785
VECTOR	pDONR223
NCBI SYMBOL	NPPB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatccccagacagacacgtccggcgctcgtctctgtcttcttgcacgtggttctgggaggtcgttccaccgcgtgggcagccccgggtcagcctcggaacttgaaacgtccgggttacaggagcagcgcaaccatttcagggcgaactgtcggagctgcaggtggagc
agacatccctggagccccccagggagagccccgtccacaggtgtctggaagtcccgagggtagccaccgagggcatccgtgggcaccgcaaatgtctctacaccctcgggcaccacgaagccccagatgggtgcaagggtctgctgttgggaggagatggaccg
gatcagctctccagtgccctgggtgcaagtgctgaggcggcattaa