

C18orf21-chromosome 18 open reading frame 21 Gene

Reference: PTXBC025950

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0711 protein C18orf21; HsT3108; PNAS-124; PNAS-131; XTP13; HBV X-transactivated gene 13 protein; HBV X-transactivated protein 13; HBV XAg-transactivated protein 13; chromosome 18 open reading frame 21		

For research use only.

GENE DESCRIPTION	chromosome 18 open reading frame 21
GENE ID	83608
ORIGIN SPECIES	Human
ACCESSIONS	BC025950
VECTOR	pDONR223
NCBI SYMBOL	C18orf21
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtgaaaaagttcaagactccaaaagtgtattgtgacactgtgaaacatgcaacagacagtgaaacatcatgtgtaaaagtagaagctttgtgcaacattgaagagcaatcctgccactcctgcaagtaaacagcctgaagacaccagagagaaggactgcaaacccaaatcat
gacatgtctggctcgaagagcaagagcccagcatcggttttcagaacacctacatctggacagtcagtagtctactgtctcctcaagaacaccagcaaaacaaagaacacttctcactaaaaatgttacttagtcagaatgaatccaaaagattccaaagggtgacttcagaaatttctt
atcttctgaagggtgacttttaaaataa