



Reference: PTXBC025989

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HLASA; LEURS; PRLTS4; mtLeuRS; leucine tRNA ligase 2, mitochondrial; leucine tRNA ligase 2, mitochondrial; leucine translase; leucyl-tRNA synthetase 2; leucyl-tRNA synthetase 2, mitochondrial		

For research use only.

GENE DESCRIPTION	leucyl-tRNA synthetase 2, mitochondrial
GENE ID	23395
ORIGIN SPECIES	Human
ACCESSIONS	BC025989
VECTOR	pDONR223
NCBI SYMBOL	LARS2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]