

LOC441046-glucuronidase, beta pseudogene Gene

Reference: PTXBC025996

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	glucuronidase, beta pseudogene
GENE ID	441046
ORIGIN SPECIES	Human
ACCESSIONS	BC025996
VECTOR	pDONR223
NCBI SYMBOL	LOC441046
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacataccggttcctccagcttcaacgacgttggccaggactggggctgcccatttttagaccagatgtggtacgaacgggaggtgaccttctggagcaatggaccaggacctgcacacaagagtggtaactgaggattgtcagtgcccactcctatgccatcgtgtgggtga
atggggctgacgcctagagcatgaggatctacctccctttgacaccgacatcagtagcctgttccagggtggggccctggccctccgcctcgcacactatcaccatcggaacatgtctatctctccacctggccaccaggagcctcctgacatggccgacacctccacgtgg
gtaccatcctgtctccaccgagacacccacctctgtgtccacccctgtggggcattacattag