

## PRRG2-proline rich Gla (G-carboxyglutamic acid) 2 Gene

Reference: PTXBC026032

|                |                                                                                                                                                                                                                                                                  |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                                                              | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                                                             | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | PRGP2; transmembrane gamma-carboxyglutamic acid protein 2; proline rich Gla (G-carboxyglutamic acid) 2; proline-rich Gla (G-carboxyglutamic acid) polypeptide 2; proline-rich Gla (G-carboxyglutamic acid) polypeptide 2; proline-rich Gla protein 2; proline-ri |                     |                |

**For research use only.**

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | proline rich Gla (G-carboxyglutamic acid) 2                                                       |
| <b>GENE ID</b>          | 5639                                                                                              |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC026032                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | PRRG2                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atgtcttgggtccagaagctcttctgggtccccagaggccagagcttctgagtagccataccggattccaagagccaaccactgggacctggagctgtcacaccagggaacctggaacgggagtgcttgaagagaggtgttctctgggaagaggccaggagattttgagga  
caacactctcacggagcgttttgggagagctacatctacaatggcaaggaggcgtggacgagtgatgtggccagcctggcgtgtgggctgacatgtggcatctgtctattgtcttggccgctgggagcctttgtgtatctgcgtggcgacagcaccgagccagcagcct  
gtcccaagaggccggctcattagccctctgagtccttgaacctcttggggccaccgacgcccctgctccacccccacccccacccaggcctccacctaagagcagcgttggcagcctctgggtacacgacgcacctccacccccctacaccagcctcaggagcctca  
ctga