

STAG3L4-stromal antigen 3-like 4 Gene

Reference: PTXBC026058

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FLJ13195; OTTHUMP00000211023; OTTHUMP00000211026; STAG3-like protein 4; stromal antigen 3-like protein 4; stromal antigen 3-like 4		

For research use only.

GENE DESCRIPTION	stromal antigen 3-like 4
GENE ID	64940
ORIGIN SPECIES	Human
ACCESSIONS	BC026058
VECTOR	pDONR223
NCBI SYMBOL	STAG3L4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagtgctgctgtagctacaagcaagaccaggatgcaggattctggagcttgttaactttttcatctgacttgtggatgtaaaggactcggggactcctctgacagctccaggtctatcctggaagaagttccagggcagcttctgtgagtttggggacattggtctgctgggtgccag
tacatcctctccatgatgacttccctatggacaacctcatctccctgctcactggtcttcagactcacaagctctgcgcttctgtcacactagcaccctggctgctatgaaactgatgacctcctggaagagttgcctccaactgagctctgcacgaagatatcaatcagctcagatgagg
ctgaaagaaacaagggccaggccaggggcacctgagcggctggagagcctgttgagaaacacaagagctgcaftaa