

ZNF585A-zinc finger protein 585A Gene

Reference: PTXBC026081

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 585A		
For research use only.			

GENE DESCRIPTION	zinc finger protein 585A
GENE ID	199704
ORIGIN SPECIES	Human
ACCESSIONS	BC026081
VECTOR	pDONR223
NCBI SYMBOL	ZNF585A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccagctaattggacctcacctcagaaatctcagccctggctccagaggatcatggcagctcctatgaggatcagtgctctcaggatgtggctatcgatttcagcagagaggatggcgccacctggacctctcagagaaacctgtaccgggatgtgatgtggagacctaca
gccacctgctctcaataggatatcaagttcctgaagcagagtggtcatgttgagcaaggaaaggaacatgggcactgcagggtgagaggccacgtcagagctgcccagcacgtgtcttgtgaactccatcacctcaagaaagcttccgagggtga