

MGC27348-ribosomal protein S2 pseudogene Gene

Reference: PTXBC026177

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	ribosomal protein S2 pseudogene
GENE ID	256355
ORIGIN SPECIES	Human
ACCESSIONS	BC026177
VECTOR	pDONR223
NCBI SYMBOL	MGC27348
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaagatcaagtcctagaggagattatctctctcgtgcccatcaaggaaatcgagatcattgactttttctggggcctctctcaaggacgaggttttgaagattatgcctgtgcggaaggagacccgcgccgccagcgccaccaggttcaaggcgtttgtgccatcagggactaca
atggctacgctggtcggttgaagtctccaaggaggtggccgcccatccgcccccatcatctgaccaagctctccattgtcaccgtgcgagaggtctactgggggaacaagatcggcaagctccacaccgtcccttgaagggtgacaggccgtgcggctctgtgctggtgc
cttcatccccgcgccagggcactggcatcgtctcagctcctgtgcccaagaagctgctcatgatagctgggtatctacgactgctacacctcagccaggggctgcactgccacctgggcaacttcgccaacgccaccttggatgctgtctctaa