



Reference: PTXBC026185

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CGI-49; NET11; saccharopine dehydrogenase-like oxidoreductase; saccharopine dehydrogenase (putative)		

For research use only.

GENE DESCRIPTION	saccharopine dehydrogenase (putative)
GENE ID	51097
ORIGIN SPECIES	Human
ACCESSIONS	BC026185
VECTOR	pDONR223
NCBI SYMBOL	SCCPDH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]