

ORM1-orosomucoid 1 Gene

Reference: PTXBC026238

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGP-A; AGP1; HEL-S-153w; ORM; alpha-1-acid glycoprotein 1; AGP 1; OMD 1; epididymis secretory sperm binding protein Li 153w; orosomucoid 1		

For research use only.

GENE DESCRIPTION	orosomucoid 1
GENE ID	5004
ORIGIN SPECIES	Human
ACCESSIONS	BC026238
VECTOR	pDONR223
NCBI SYMBOL	ORM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtgtcctgggtttctacagctcctcctacctcctggaagccagatccattgtgtgccaacatgacccgtgcccacccaacgccacccctggaccggatcactggcaagtgttttatatcgcacgtggttgaacacgaggtacataaagtcggtcaggaga
tccaagcaaccttcttttacttcccccaacaagacagaggacagatctttctcagagagtaccagaccgacagggaccagtgcactataaacaccacctacctgaatgtccagcgggaaafgggaccatctccagatagctggagggccaagagcatttcctcacttgctgacctca
gggacaccaagacatcatgcttctttgacgtgaacgatgagaagaactggggcctgtctgtctatgctgacaagccagagacaccaaggagcaactgggagagttctacgaagctctcactgcttgcgcattcccaagtcagatgtctgtacaccgattgaaaaaggataagtg
tgagccactggagaagcagcagcagaaggagaggaacaggaggagggggaatcctag