

MAPKSP1-MAPK scaffold protein 1 Gene

Reference: PTXBC026245

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MAPKSP1; MAP2K1IP1; MAPBP; MP1; PRO0633; Ragulator3; ragulator complex protein LAMTOR3; MAPK scaffold protein 1; MEK binding partner 1; late endosomal/lysosomal adaptor and MAPK and MTOR activator 3; mitogen-activated protein kinase kinase 1 interacting p		

For research use only.

GENE DESCRIPTION	MAPK scaffold protein 1
GENE ID	8649
ORIGIN SPECIES	Human
ACCESSIONS	BC026245
VECTOR	pDONR223
NCBI SYMBOL	MAPKSP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggatgacctaagcgattctgtataaaaagttaccaagtgttgaagggtccatgccattgtgtgcagatagagatggagtacctgttataaaggaggcaaatgacaatgctccagagcatgtttacgacctggttcttaccactttgcccttgcaacagaccaaggaagcaaac
ttggactttccaaaaataaaagtatcatctgttactataacacctaccaggtgttcaatttaacgtttacctttgtgtgtgagtttcatagccagcagcagtgccaatacaggactaattgtcagcctagaaaaggaaactgtccactgtttgaagaactgagacaagttgtggaagtctttaa