

INSL4-insulin-like 4 (placenta) Gene

Reference: PTXBC026254

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	EPIL; PLACENTIN; early placenta insulin-like peptide; early placenta insulin-like peptide (EPIL); insulin-like 4 (placenta); insulin-like peptide 4; insulin like 4		

For research use only.

GENE DESCRIPTION	insulin-like 4 (placenta)
GENE ID	3641
ORIGIN SPECIES	Human
ACCESSIONS	BC026254
VECTOR	pDONR223
NCBI SYMBOL	INSL4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccagcctgttccggctctatctgccagcaatctggctgctgtgagccaactccttagagaagcctagcagcagctgaggggatgtgtcccgattggaaaacactgtgtcatattgcccatgctgagaagacattcaccacccaggggtgctgtggaatc
tggacgtcccaagaagaatggtgtcaactccaacaagaagatggacaagccttaggtacgacatcagaattccttaattgtcaccagagctgaagaaccactgtctgaaggcagccatcattgaagaaaataatacttcccgcaaaaagagaagtggacgtcacagattgatcc
attctgtgtgaagtaattgtgacgatggaactcagtttaattatgtacatag