

C7orf11-chromosome 7 open reading frame 11 Gene

Reference: PTXBC026265

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C7orf11; ABHS; ORF20; TTD4; M-phase-specific PLK1-interacting protein; Russell-Silver syndrome region; TTD non-photosensitive 1 protein; M-phase specific PLK1 interacting protein		

For research use only.

GENE DESCRIPTION	chromosome 7 open reading frame 11
GENE ID	136647
ORIGIN SPECIES	Human
ACCESSIONS	BC026265
VECTOR	pDONR223
NCBI SYMBOL	C7orf11
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcagcgcagagaattttcgcccccaactctcttaccctgtccgggtgaggaggttgggtagcggagcagcttccgggaaccccgggcgggcggggaccacggccgctccctcgagacgggtacgggagtcgcaccacacgcgcgtacgggccccgtcta
ggcgtacgggagcagtcactctccgcgacagggcgagcttccggggggcggttcgggtctccgtccctggcgctacctggctctactccaggtcccccggggtcccgagcagcaattcgctactccccagggcagcagcagacccccaggggtctccaagg
acatctacacatttggatcaggcggtgttagagaaaaagaatgtctaagtgtgaaaattattcaagccttcaatgctgaagatcctgggctggcctagaaccagtatctgtatggatataagccaacaatacagaatactcaaacattcacaggcaaaaagggaagatactttgtt
aa