

## BZW1-basic leucine zipper and W2 domains 1 Gene

Reference: PTXBC026303

|                |                                                                                                                                                |                     |                |
|----------------|------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                            | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                           | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | Nb1a10236; basic leucine zipper and W2 domain-containing protein 1; basic leucine-zipper protein BZAP45; basic leucine zipper and W2 domains 1 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | basic leucine zipper and W2 domains 1                                                            |
| <b>GENE ID</b>          | 9689                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC026303                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | BZW1                                                                                             |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atgaataacaaagcagcaaaagccaacgctatcaggccagcgttttaaaactagaaaaagagatgaaaaagagaggtttgacccactcagttcaagactgtattattcaagccttaactgaaacccggtactgatttgaagcagtagtaagttcttgatgcttctggagcaaaactgat  
taccgtcgatgacagaacactctttgacattctggtggctgggtgaatgctggccccaggtgtacactggcagatgacatgatgcgtacagatgtctgctgtttgagcccaagaagatctagagaccatgcaagcatttgcacaggttttaacaagttatcaggcgtacaaatcct  
ggagaaaggttttgaagatgaagtaaaaaagctgctgctgttcctgaagggttttcagagctggagaggaacaagctagctatgttgactggtgttcttctggcctaaggaacacttaagcatccattcttaatagcctttataatgaaaatttggtaagaagaggttcagcagcttttgcgtg  
aagctctttaaatcatggataaatgaaaaagataatcaatgcagtagctgcaagctctgggaaagtcagcatggataacagactgatggaaactcttctgccaataagcaaaagtgttgaacacttcacaaaatatttactgaggcaggcttgaagagcttcagaatatttcgggaatcagcaa  
accatcggagctgtaaggagctccagaagaacttcaagaacagatgcccgtggatccatttaaggatataattttatatgtcaaggagagatgaaaaaaacaacatcccagagccagttgtcatcggaaatagctgtgcaagtgtaagcactgtggaatggaacaaaaaga  
ggagctttagcagagcaagccatcaagcactgaagcaatacagccctacttgcctttactactcaaggctcagctgagctgactctgttactgaaattcaggagattgtatgacaacattcattcatgaagcctccagaaaatagtggtgctttttataaagctgaagtcctga  
gcgaggagccatttgaagtggtataaagatgcacatgttgcagggggaagagtggttcttgagcaaatgaaaaagttttagaatggctcaaaaatgctgaagaagaatctgaatctgaagctgaagaaggtgactga