



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SLC48A1-solute carrier family 48 (heme transporter), member 1 Gene

Reference: PTXBC026344

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

HRG-1; HRG1; hHRG-1; heme transporter HRG1; heme-responsive gene 1 protein homolog; solute carrier family 48 (heme transporter), member 1; solute carrier family 48 member 1

For research use only.

GENE DESCRIPTION	solute carrier family 48 (heme transporter), member 1
GENE ID	55652
ORIGIN SPECIES	Human
ACCESSIONS	BC026344
VECTOR	pDONR223
NCBI SYMBOL	SLC48A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcccccgtcagggcgtgcagctcggccctccgcgcgcctactcgcgcacagctcgcggcgccctctccactctctcgtctgcacgggtgtctaccgacagccggggaccggcccatggaggggctcgcaggcgggtctgcgacgtgtggctctctgtagcgcacgtgatgtacat
ggcagatatttggaggacgtgcctcaaggggctgcgcgcgtctctctctggcgcgtctctctcgcgcgcctctcgcctcgtcgcctctcgcacctctcgtctgcacctcctcgtctgcacacaccgggcacagcctcacagaccgccacagctactactctcagcgtctggagcctcaattctctcaagtg