

## KIF19-kinesin family member 19 Gene

Reference: PTXBC026362

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	kinesin-like protein KIF19; KIF19A; kinesin family member 19		

**For research use only.**

---

<b>GENE DESCRIPTION</b>	kinesin family member 19
<b>GENE ID</b>	124602
<b>ORIGIN SPECIES</b>	Human
<b>ACCESSIONS</b>	BC026362
<b>VECTOR</b>	pDONR223
<b>NCBI SYMBOL</b>	KIF19
<b>LENGTH</b>	ProteoGenix bp
<b>SEQUENCE PRIMERS</b>	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

### ORF (Only for reference, product provided without STOP codon)

atgagcagatggtgttctcatggaccaatggagatcccacgacatcctgcggcgcatcgtctccgggagagtcctacctgttcgacgtggcctttgacttcaccgccaccaggagatggtgtatcaggccaccaccaagacctcatcgaggcgctcatctcaggtacaatg  
ccactgtctttgcctatgccccacagctgtgggaaaaacctacacctgctgggcacagaccaggagcctggcctctatgttcagaccctcaacgacctctccgtgccatcgaggagaccagcaatgacatggagtatgaggtctccatgtctactcctggagatcatgcagctgctga