

PRPF40A-PRP40 pre-mRNA processing factor 40 homolog A (S. cerevisiae) Gene

Reference: PTXBC027178

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FBP-11; FBP11; FLAF1; FBNP3; HIP-10; HIP10; HYPA; NY-REN-6; Prp40; pre-mRNA-processing factor 40 homolog A; Fas-ligand associated factor 1; Huntingtin-interacting protein A; NY-REN-6 antigen; PRP40 pre-mRNA processing factor 40 homolog A; fas ligand-assoc		

For research use only.

GENE DESCRIPTION	PRP40 pre-mRNA processing factor 40 homolog A (S. cerevisiae)
GENE ID	55660
ORIGIN SPECIES	Human
ACCESSIONS	BC027178
VECTOR	pDONR223
NCBI SYMBOL	PRPF40A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaaacgaaagaatctgcatltaagatgtttaaaacaagctgctcctccgataaattggatgctgtctgggaagatccgtgagagatttgaagagccagcatttggagacataactctagaatctgaagaaaaacgaatatttaagattttatgcatgtgcttgagcatgaatgtagcatcatcattcaaaagaacaagaacattctaaagaatctaaaaaacatcataggaacgttcccgtctcgatcggttcagattcagatgatgatagccattcaaaagaaaaaacagcagcagatcagagtcgtctgtctcagaaacattctctagtgcagagctgagagaagttataaaaagtcaaaaagcataagaagaaaaagtaagaagagagacataaatctgactctccagaatccgatgctgagcgcagagaaaggataaaaaagaaaaagatcgggaaagtgaagaaagacagaactagacaaagatcagaatcaaaacacaaatcgcttaagaaaaagactggaaaaggattctgtaattgggatactctggcagcgaactgagtgaaggggaattgaaaagccagagaagaacctttggagcaactggatgatcatcaataa