



Reference: PTXBC027241

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MER57B-ALDH1L1; 10-FTHFDH; 10-fTHF; FDH; FTHFD; cytosolic 10-formyltetrahydrofolate dehydrogenase; 10-formyltetrahydrofolate dehydrogenase; aldehyde dehydrogenase 1 family member L1		

For research use only.

GENE DESCRIPTION	aldehyde dehydrogenase 1 family, member L1
GENE ID	10840
ORIGIN SPECIES	Human
ACCESSIONS	BC027241
VECTOR	pDONR223
NCBI SYMBOL	ALDH1L1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr