



Reference: PTXBC027332

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	membrane protein FAM174A; 2310044D20Rik; Fam174; Naa6; Tmem157; transmembrane protein 157; family with sequence similarity 174, member A		

For research use only.

GENE DESCRIPTION	family with sequence similarity 174, member A
GENE ID	345757
ORIGIN SPECIES	Human
ACCESSIONS	BC027332
VECTOR	pDONR223
NCBI SYMBOL	FAM174A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

atgaagcctgcctgactgctgctgctgctcagccacctcttgcttcctgctcctctctgctgttgctgctgctgaactaagccggccctcgagctcgtgctgcagcagccgaggccgccaggctctgggctcctgacctagaccacgacattaccgcctgccaccggcccttaccctgccccagcagccggccgtgctgtgctggaagctgcggggccgggggctcgcaggaggcaalggcgagcaacctgtggccgggcttgagacggagcatcaggaggggaaagccgggggaagcctcggtgggtggccgctgctgtgagccccaacctgtgcgacaaagccatgaccagccggcgctgacagctggtggtgagcggcgcggtgctgctggtgactctgcgcagcagcgtgcaggtgagaaagacccaagaagactaggagaatagtgagtttggacactaacatagaaaatgtgaattgacaccttttagaacaggaatgatgaggatga